


How an MCP works



## Nickel-chromium electrode

For over thirty years, PHOTONIS (formerly Philips Photonics, Galileo and Burle) has consistently set the standard in electron multipliers and related products. Today, an extensive R&D program coupled with unsurpassed expertise in microchannel plate (MCP) technology continues to deliver a succession of product and process improvements that push aside previous technology limits. This expertise is used to mass-produce MCPs for our image intensifier tubes as well as for a variety of custom scientific applications. This, plus the ability to meet your product requirements, makes Photonis the preferred choice of professionals the world over.

A microchannel plate (MCP) is an electron multiplier for detecting X-rays, ultraviolet radiation and charged particles. The output is a two-dimensional electron image which preserves the spatial resolution of the original input radiation, but with a linear gain up to 1000. This may be used for exciting a phosphor screen placed close to the output, giving a visual representation of the radiation pattern. Alternatively, the electron image can be read out by, for example, a wedge-and-strip or fast delay-line anode array.

Important features of MCPs are:

- high electronic gain
- immunity from magnetic fields
- fast response
- low noise
- low power consumption
- high spatial resolution
- small size and ruggedness

Each plate consists of an array of tiny glass tubes fused together to form a thin disc. Both faces of the disc are metal-coated to provide parallel electrical connections to all channels. In a vacuum, and with a potential difference (usually 800 to 1400 V) across the plate, each channel becomes a continuous dynode electron multiplier, operating on the same principle (electron avalanche) as its cousin - the single-channel electron multiplier.

MCP

Special MCP's

Micropore optics

Photonis has developed the first-ever square-pore radially-packed X-ray focusing MCP optics and is currently refining the technology with continued support by ESA. For special projects, options are available, such as a square-pore square-packed X-ray optics with pore sizes of 10 to 100  $\mu\text{m}$ , plates with a thickness of several millimeters, formed with a spherical radius or with metal-coated pore walls for optimum reflectivity can be developed and supplied.

Curved-channel MCPs

The curved-channel MCP was invented by the Philips LEP laboratory 1973 and produced in our factory in Brive in 25  $\mu\text{m}$  pore, 25 mm diameter for some time. It was seen as a new way to prevent ion feedback at high gain with preserved spatial resolution. The major application was invention of the Multi-Anode Multi-Array (MAMA) tube by Timothy et al. at the time. Today, improved MCP readout arrangements in the form of "wedge & strip" anodes, fast delay-lines or Vernier anodes offer equal or better spatial resolution with the MCP stacks mentioned above. As the curved-channel MCPs were very difficult and costly to make in large formats they are only used for a few MAMA tubes, offered by competitors who copied the process after our patents expired.